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
1. Furukawa M, Minegishi Y, Izumiyama S, Yagita K, Mori H, Uemura T, Etoh Y, Maeda E, Sasaki M, Ichinose K, Harada S, Kamata Y, Otagiri M, Sugita-Konishi Y, Ohnishi T: The development of a novel, validated, rapid and simple suitable method for the detection of *Sarcocystis fayeri* in horse meat in the sanitary control setting. *Biocontrol Science*, **21**, 131-134 (2016)

EasyAmp *Sarcocystis fayeri* Detection Kit

Code No. 391-5310

I Product Description

This product is a DNA test kit to detect *Sarcocystis fayeri* by using Loop-mediated Isothermal Amplification (LAMP) method. 2×LAMP Master Mix contains isothermally amplifiable strand displacement DNA polymerase, inorganic pyrophosphatase, Mg²⁺, dNTPs, optimum buffer, and fluorescent dye that binds to the double-stranded DNA and allows detection of amplified DNA by fluorescence device.

[Precaution]

1. This product is a kit to detect *Sarcocystis fayeri* and cannot be used for other purpose.
2. This product should be handled only by persons having a basic knowledge about reagents.
3. This product should be handled in accordance with the descriptions in the manual. We are not responsible for problems caused if this product is not handled in accordance with the manual.
4. Expire date is printed on the label.
5. Follow the instructions from your region to discard this product.
6. Do not electrophorese and autoclave amplification product by LAMP to prevent contamination of test environment.
7. Read SDS before use.
SDS can be downloaded from our web site (<https://nippongene.com/kensa/list.html>).
8. LAMP (Loop-mediated Isothermal Amplification) is the isothermal gene amplification technique developed by Eiken Chemical Co., Ltd.
9. Turbidity measuring Real-time device for LAMP (Eiken Chemical Company) cannot be used for this kit.

II Contents

For 50 reactions

	Contents	Volume	Qty
I	2×LAMP Master Mix	640 μL	1
II	<i>Sarcocystis</i> Positive Control	150 μL	1
III	<i>Sarcocystis</i> Primer Mix	375 μL	1
IV	Manual	—	1

Storage

- ◆ Store at -20°C and dark place.
Use before expire date.
Residual reagent should be stored at -20°C again.

III Necessities in addition to this product

- Micropipette and tip
- Centrifuge
- Micro tube and rack
- Heat block
- DNA extraction reagent
Recommend: EasyPrep DNA Extraction Reagent (NIPPON GENE, Code No. 391-4310)
- fluorescence detecting device
example) **Genie**[®] II, **Genie**[®] III, Real-time PCR device
- Pestle
- Vortex
- Crush ice
- Sterilized water

IV Protocol

1. Extraction of DNA from horse muscle

important To prevent the contamination of test environment, please extract DNA in isolated place from the place this product is used.

1-1. Protocol using EasyPrep DNA Extraction Reagent (separately sold)

- 1) Weigh 300 mg of horse muscle to 1.5 mL tube
- 2) Add 300 μ L of EasyPrep DNA Extraction Reagent and grind and mix with pestle, and then vortex well.
- 3) Incubate for 10 min at 95°C.
- 4) Incubate for 1 min on ice.
- 5) Centrifuge at 15,000xg for 5 min at room temperature.
- 6) To prepare template DNA, transfer supernatant into a new 1.5 mL tube and dilute to 5 times by sterilized water.

example: Transfer 100 μ L of supernatant into a new 1.5 mL tube contained 400 μ L of sterilized water.

2. Preparation of LAMP reaction solution

2-1. Preparation of reagents

Melt the reagents completely at room temperature. Vortex the reagents to dissolve uniformly and then spin down and leave them to stand on ice.

2-2. Preparation of LAMP reaction solution

- 1) Prepare Pre Mix .

Reagents	1 reaction	8+1 reactions
2xLAMP Master Mix	12.5 μ L	112.5 μ L
Sarco Primer Mix	7.5 μ L	67.5 μ L
Total	20.0 μ L	180.0 μ L

- 2) Dispense the 20 μ L of Pre Mix.
- 3) Add template DNA to the Pre Mix.

Add 5 μ L of negative control (sterilized water) and close the cap of the tube. Same as negative control add 5 μ L of template DNA (see 1-1) and close the cap of the tube. At last, add 5 μ L of positive control and close the cap of the tube.

3. LAMP reaction and measuring

Mix the LAMP reaction solution by inverting and spin down. Set to the fluorescence detecting device. Start measuring (DNA amplification monitoring and Annealing curve analysis or Melting curve analysis). Follow the setting method of the device's manual.

3-1. LAMP reaction

66°C for 30 min → Annealing curve analysis or Melting curve analysis

3-2. Measuring method

Set Fluorescent wavelength to ResoLight Dye or SYBR™ Green I when using Real-time PCR device.

4. Determination result

Test for detection of *Sarcocystis fayeri* is accepted when all of the following requirements are satisfied.

- Amplification of negative control is not detected.
- Amplification of positive control is detected.

The result is determined positive when all of the following requirements are satisfied.

- Amplification of template DNA is detected.
- In Annealing curve analysis or Melting curve analysis, the temperature of template DNA is within ± 1 °C of the temperature of positive control.

V Reference

1. Furukawa M, Minegishi Y, Izumiyama S, Yagita K, Mori H, Uemura T, Etoh Y, Maeda E, Sasaki M, Ichinose K, Harada S, Kamata Y, Otagiri M, Sugita-Konishi Y, Ohnishi T: The development of a novel, validated, rapid and simple suitable method for the detection of *Sarcocystis fayeri* in horse meat in the sanitary control setting. *Biocontrol Science*, **21**, 131-134 (2016)