GM Soy $\Delta\Delta$ Cq Standard Plasmid Set

Code No. 316-09131 (1 Set)

Storage:

Store at -20°C, or at -80°C for long term storage.

Component:

This kit contains the standard plasmid DNA used for qualitative testing by the $\Delta\Delta$ Cq method using real-time PCR for soybean grains, and blank sample (ColE1 DNA/TE) as non-template control.

ColE1 DNA/TE	(ColE1 DNA 5 ng/ μ l)	180 μ l \times 2
GM Soy $\Delta\Delta$ Cq Standard Plasmid 50	(50 copies/ μ l)	120 μ l \times 2
GM Soy $\Delta\Delta$ Cq Standard Plasmid 100K	(100,000 copies/ μ l)	60 μ l \times 2

Usage example:

Reaction component (Total volume 25 μ l/well)	
TaqMan [®] Universal Master Mix (2X)	12.5 μ l
Primer pair (25 μ M each of F&R)	0.5 μ l
Probe (10 μ M)	0.5 μ l
Standard Plasmid DNA	5 μ l
Water	6.5 μ l

Note:

One package of the kit is controlled the quality as one product. We cannot guarantee qualitative test result by the $\Delta\Delta$ Cq method if reagents from other packages are used together.

Storage precautions:


- Low concentration DNA is unstable. We recommend storing the kit at -80°C for long-term storage.
- Avoid repeated freeze-thaw cycles more than 6 times.
- Do not aliquot the plasmid for storage to avoid a risk of contamination.
- All reagents should be stored under the same conditions.

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.
2. Takabatake R, Onishi M, Koiwa T, Futo S, Minegishi Y, Akiyama H, Teshima R, Furui S, Kitta K.: *Shokuhin Eiseigaku Zasshi*, **51** (5), 242 (2010)

License:

This product is sold under licensing arrangements with National Food Research Institute, Asahi Breweries, Ltd., Nippon Flour Mills Co., Ltd.

Primer pair and Probe name:

- RR2** GM Soybean (RR2) Detection RR2 Oligonucleotide Set (Code No. 319-07661)
- Le1** GM Soybean (RRS) Detection Le1 (for endogenous gene) Oligonucleotide Set (Code No. 319-05601)
- P35S** GMO Detection P35S Oligonucleotide Set (Code No. 313-05621)

Figure: Structure of GM Soy $\Delta\Delta$ Cq Standard Plasmid (pSCS)

This product is a linearized plasmid diluted with ColE1 DNA/TE Solution.

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