# GeneAce RL qPCR Mix with UNG

## I. Description

GeneAce RL qPCR Mix with UNG is a real time PCR reagent for the intercalator method. The Dye contained in this product generates a higher fluorescence signal compared to other fluorescence dyes and PCR products can be analyzed with higher resolution in melting curves.

Because this product contains dUTP instead of dTTP, Carryover contamination can be prevented by using included UNG. UNG is an enzyme that degrades PCR products containing uracil bases by reacting at 50°C for 2 minutes before the PCR reaction.

50x ROX Passive Reference might be required for using real time PCR instruments that normalize the fluorescent signal by ROX (This product does not contain ROX).

## II. Storage

-20°C, protected from light.

This product should be protected from light. After thawing, the master mix may be stored at 4°C for up to one month or returned to -20°C for long term storage by expired date.

## III. Components

**200 rxns**

2 x GeneAce RL qPCR Mix No ROX *1)	1250 µL x 4
------------------------------------	-------------

Uracil-DNA Glycosylase (UNG), 1U/ µL	100 µL x 1
--------------------------------------	------------

ddWater	1250 µL x 4
---------	-------------

\*1) It contains Hot-Start Gene Taq NT, dNTP Mixture with dUTP, Mg<sup>2+</sup>, fluorescence dye, and stabilizer

## IV. Notes

- **The PCR must start with an initial 10 minutes incubation at 95°C to activate the chemically modified Hot-Start Taq DNA polymerase.**
- **We recommend two-step protocols. (See V. Typical PCR protocol)**
- Always ensure that the product has been fully thawed and mixed before use. Gently mix the mixtures without creating bubbles.

## V. Typical PCR Protocol

Choose an appropriate total reaction volume, depending on the instrument used. (e.g., 50 µL or 25 µL)

Choose ResoLight Dye or SYBR® Green I for the fluorescence wavelength setting.

The Dye contained in 2x GeneAce RL qPCR Mix is excited by light with a wavelength between 450 nm and 500 nm and the maximum excitation is at 487 nm. The fluorescence wavelength is between 480 nm and 560 nm and the maximum is at 503 nm.

Component	Volume	Volume
2 x GeneAce RL qPCR Mix No ROX	25.0 µL	12.5 µL
25 µM each Primers	1.0 µL	0.5 µL
Uracil-DNA Glycosylase (UNG)	0.5 µL	0.25 µL
Template	5.0 µL	2.5 µL
d.d.H <sub>2</sub> O	up to	50.0 µL 25.0 µL

Recommended cyler conditions

<b>50°C</b>	<b>2 min.</b>	<b>UNG treatment</b>
<b>95°C</b>	<b>10 min.</b>	<b>Enzyme activation</b>
<b>95°C</b>	<b>30 sec.</b>	) <b>45 cycles</b>
<b>60°C</b>	<b>1 min.</b>	

**Melting curve analysis**

Optimum PCR condition depends on primer sequence and template DNA and so on.

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