GeneAce SYBR™ One Step RT-qPCR Mix

Code No. 319-09741, 315-09743

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Store at -20°C(遮光)

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Store at -20°C (protected from light)

Description:

GeneAce SYBR™ One Step RT-qPCR Mix is a convenient two-component reagent designed for real-time RT-qPCR using SYBR™ Green I dye.

This product ensures high specificity, reproducibility and high amplification efficiency through optimized buffer components, Hot-Start *Taq* DNA Polymerase, and a modified M-MLV Reverse Transcriptase. It also includes a unique passive reference dye compatible with a wide range of instrument platforms.

Components:

(1) 2X GeneAce SYBR™ One Step Mix

1.25 mL x 1 (125 rxns)

1.25 mL x 4 (500 rxns)

Contains: Hot-Start *Taq* DNA Polymerase,
dNTP, Mg²⁺, SYBR™ Green I,
Passive Reference Dye,
Stabilizers

(2) 40X RT Mix (for SYBR™)

63 µL x 1 (125 rxns)

63 µL x 4 (500 rxns)

Contains: Modified M-MLV Reverse
Transcriptase, RNase Inhibitor

Notes:

- Due to the use of an antibody-bound hot-start *Taq* DNA polymerase, an initial denaturation step at 95°C for 2 minutes is required to activate the enzyme.
- A two-step cycling protocol is recommended. (See Typical PCR protocol.)
- Always ensure that the product is completely thawed and thoroughly mixed before use. Mix gently to avoid generating bubbles. During use, keep this component on ice. After use, store immediately at -20°C.
- Only specific primers are available for reverse transcription reaction of this product. Do not use random primers or oligo(dT) primers.
- Uracil-N-Glycosylase treatment is not applicable to this product.

Typical PCR Protocol:

Choose an appropriate total reaction volume, depending on the instrument used.

Component	Volume [20 µL rxn]
2x GeneAce SYBR™ One Step Mix	10.0 µL
40x RT Mix (for SYBR™)	0.5 µL
25 µM each Primers	0.4 µL
Template	2.0 µL
d.d.H ₂ O	up to 20.0 µL

Recommended thermal cycling conditions. *

50°C 10 min. Reverse transcription

95°C 2 min. Enzyme activation

95°C 10 sec.
60°C 45 sec.) 45 cycles

Melting curve analysis

- * Optimal PCR condition depends on primer sequence and template RNA and so on.
- * When using Thermo Fisher Scientific real-time PCR instruments, you can select either the "Standard" or "Fast" run mode. Choose the run mode based on your experimental purpose. Be sure to change the setting from the default to the "Recommended thermal cycling conditions," and ensure that both the reverse transcription step and the enzyme activation step are included.
- * If the amplicon length is 200 bp or longer, it is recommended to increase the extension time from 45 to 60 seconds. Design primers so that the amplicon length does not exceed 300 bp. Ideal length is 80–150 bp.

For Licensee Products sold in the Research Field:

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.