



GM Maize Detection MIR604 Oligonucleotide Set

Code No. 311-08341 (1 Set)

Storage:

Store at -20°C and protect from light.

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

Detection of:

MIR604

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.
2. Mano J, Furui S, Takashima K, Koiwa T, Futo S, Minegishi Y, Akiyama H, Teshima R, Kurashima T, Takabatake R, Kitta K.: *Food Hygiene and Safety Science (Shokuhin Eiseigaku Zasshi)*, **53** (4), 166 (2012)

Probe & primer name:

MIR604 primer F & MIR604 primer R

MIR604 probe (FAM/TAMRA)



GMトウモロコシ系統別 DNA MIR604 オリゴヌクレオチドセット

Code No. 311-08341 (1 Set)

保存:

Store at -20°C(遮光)

内容:

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

検出:

MIR604

プローブ&プライマー名:

MIR604 primer F & MIR604 primer R

MIR604 probe (FAM/TAMRA)

参考文献:

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」
2. Mano J, Furui S, Takashima K, Koiwa T, Futo S, Minegishi Y, Akiyama H, Teshima R, Kurashima T, Takabatake R, Kitta K.: *Food Hygiene and Safety Science (Shokuhin Eiseigaku Zasshi)*, **53** (4), 166 (2012)

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。