



GMO Detection P35S Oligonucleotide Set

Code No. 313-05621 (1 Set)

Storage:

Store at -20°C

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 90 µL

Detection of:

P35S

Probe & primer name:

P35S-1 [P35S 1-5' & P35S 1-3']

P35S-Taq (FAM/TAMRA)

Amplicon length:

101 bp

Reference:

1. "JAS Analytical Handbook, Manual of Inspection and Analysis for GM-Foods"; in Japanese, by IAI Center for Food Quality, Labeling, and Consumer Services of Japan.

2. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.

License:

This product is sold under licensing arrangements with National Food Research Institute, Asahi Breweries, Ltd., Nippon Flour Mills Co., Ltd.

組換え DNA P35S オリゴヌクレオチドセット

Code No. 313-05621 (1 Set)

保存:

Store at -20°C

内容:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 90 µL

検出:

P35S

プローブ&プライマー名:

P35S-1 [P35S 1-5' & P35S 1-3']

P35S-Taq (FAM/TAMRA)

増幅長:

101 bp

参考文献:

1. 農林水産消費技術センター「JAS 分析試験ハンドブック 遺伝子組換え食品検査・分析マニュアル」
2. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」

License:

本品は国立研究開発法人 農業・食品産業技術総合研究機構、アサヒビール株式会社、日本製粉株式会社により、特許第 4291568 号として登録されています。本品を有償での受託検査等の営利目的に用いる場合は、別途ライセンスが必要になります。

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。