



GMO Detection NOS Oligonucleotide Set 2

Code No. 310-08551 (1 Set)

Storage:

Store at -20°C

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

Detection of:

TNOS (NOS terminator)

Probe & primer name:

NOS ter 3-5' & NOS ter 2-3'

NOS-Taq (FAM/TAMRA)

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.

License:

This product is sold under licensing arrangements with National Food Research Institute, Asahi Breweries, Ltd., Nippon Flour Mills Co., Ltd.

組換え DNA NOS オリゴヌクレオチドセット 2

Code No. 310-08551 (1 Set)

保存:

Store at -20°C

内容:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

検出:

TNOS (NOS terminator)

プローブ&プライマー名:

NOS ter 3-5' & NOS ter 2-3'

NOS-Taq (FAM/TAMRA)

参考文献:

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」

License:

本品は国立研究開発法人 農業・食品産業技術総合研究機構、アサヒビール株式会社、日本製粉株式会社により、特許第 4291568 号として登録されています。本品を有償での受託検査等の営利目的に用いる場合は、別途ライセンスが必要になります。

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。