



GMO Detection PAT Oligonucleotide Set

Code No. 317-09781 (1 Set)

Storage:

Store at -20°C and protect from light.

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

Detection of:

PAT

Probe & primer name:

PAT-120F & PAT-120R

PAT-1-Taq (FAM/TAMRA)

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.

組換え DNA PAT オリゴヌクレオチドセット

Code No. 317-09781 (1 Set)

保存:

Store at -20°C (遮光)

内容:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

検出:

PAT

プローブ&プライマー名:

PAT-120F & PAT-120R

PAT-1-Taq (FAM/TAMRA)

参考文献:

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。