



# **GM Papaya Detection PRSV-YK(YK-2) Oligonucleotide Set**

---

**Code No. 313-81201** (1 Set)

**Storage:**

Store at -20°C

**Component:**

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

**Detection of:**

PRSV-YK

**Probe & primer name:**

YK-2F & YK-2R

YK-2P (FAM/TAMRA)

**Usage example:**

Use 0.8 µL of this product (Primer pair) in each 25 µL reaction.

For detail on the PCR conditions, refer to the References.

**Reference:**

1. "Testing for Foods Produced by Recombinant DNA Techniques" by Ministry of Health, Labour and Welfare, Japan.

# GM パパイヤ系統別 DNA PRSV-YK(YK-2) オリゴヌクレオチドセット

Code No. 313-81201 (1 Set)

**保存:**

Store at -20°C

**内容:**

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

**検出:**

遺伝子組換えパパイヤ 1 系統 (PRSV-YK)

**プローブ&プライマー名:**

YK-2F & YK-2R

YK-2P (FAM/TAMRA)

**使用例:**

遺伝子組換えパパイヤ系統 (PRSV-YK) 検知試験プライマー対として PCR 系 25 μL に本品 (Primer pair) を 0.8 μL 使用する。

PCR 条件の詳細については、参考文献を参照してください。

**参考文献:**

1. 厚生労働省「安全性未審査の組換え DNA 技術応用食品の検査方法」

本品は、試薬 (試験研究用) として販売しているものです。  
医薬品の用途には使用しないでください。