



GM Papaya Detection 55-1 Positive Control Plasmid 2

Code No. 317-07721 (50 μ L)

Code No. 313-07723 (50 μ L \times 5)

Storage:

Store at -20°C

Component:

Plasmid 50 μ L

Usage example:

Use 5 μ L of this plasmid in each 25 μ L reaction (e.g. Add 5 μ L of this plasmid to 20 μ L of PCR mixture.)

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.




Figure: Structure of GM Papaya Detection 55-1 Positive Control Plasmid 2

This product is a linearized plasmid diluted with ColE1 DNA/TE Solution.

GM パパイヤ系統別 DNA 55-1 陽性コントロールプラスミド 2

Code No. 317-07721 (50 μ L)

Code No. 313-07723 (50 μ L \times 5)

保存:

Store at -20°C

内容:

Plasmid 50 μ L

使用例:

PCR 系 25 μ L の場合、本品を 5 μ L 使用する(e.g. PCR 混合液 20 μ L に本品を 5 μ L 添加する。)

参考文献:

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。




Figure: GM トウモロコシ系統別 DNA Bt10 陽性コントロールプラスミドの構造

本品は、ColE1/TE 溶液によって希釈された直鎖状プラスミド DNA です。